



# AlerTox<sup>®</sup> Sticks

Rapid tests for allergen detection in food, drinks, and surfaces for internal quality control.



## AlerTox<sup>®</sup> Sticks advantages:

- Quick and Reliable
- Results in 10 minutes
- No Cross-Reactivity
- High Sensitivity

# AlerTox<sup>®</sup> Sticks are a practical tool for the control of allergen content in raw materials, final products, and working surfaces.

## Industrial application

- Semi-quantification of allergens in food
- Food quality control
- Trace allergen contamination in food, liquid and surfaces
- Implementation and fulfillment of HACCP, ISO 22000, BRC, and IFS compliance



Product Name	Ordering No.	Quantity
AlerTox <sup>®</sup> Sticks Egg	KIT3025	10 tests
AlerTox <sup>®</sup> Sticks $\beta$ -Lactoglobulin	KIT3018	10 tests
AlerTox <sup>®</sup> Sticks Casein	KIT3021	10 tests
AlerTox <sup>®</sup> Sticks Total Milk	KIT3015	10 tests
AlerTox <sup>®</sup> Sticks Crustacean	KIT3036	10 tests
AlerTox <sup>®</sup> Sticks Fish	KIT3038	10 tests
AlerTox <sup>®</sup> Sticks Mustard Seed	KIT3040	10 tests
AlerTox <sup>®</sup> Sticks Soy	KIT3028	10 tests
AlerTox <sup>®</sup> Sticks Hazelnut	KIT3035	10 tests
AlerTox <sup>®</sup> Sticks Almond	KIT3033	10 tests
AlerTox <sup>®</sup> Sticks Walnut	KIT3092	10 tests
AlerTox <sup>®</sup> Sticks Peanut	KIT3094	10 tests
AlerTox <sup>®</sup> Sticks Sesame	KIT3096	10 tests

**NOTE:** As with any qualitative or semiquantitative test, results should be periodically validated with a quantitative test performed in an outside laboratory.